

Passive photodiode visibility — run report

RUN

214651_pd-visibility_pd-visibility

Status: failed | Started: 2026-09-21T20:46:51.133184+00:00

PURPOSE

Measure fringe contrast on the selected PD input under external drive. No RP output initialization or PAX connection.

REQUIRED OPTICAL SETUP

DC-coupled PD on the selected RP input; external drive must span complete fringes. Set its frequency and, when known, the blocked-light voltage offset. Existing Pyrpl monitor server required.

RECORDED BENCH CONTEXT (RECIPE)

In front of NPBS 1, C port

Red Pitaya (0-1 V command) -> Thorlabs MDT690 -> piezo

DATA

0 recorded samples

Source: data.csv; full acquisition settings: recipe.json; status: run.json

RUN ISSUE

PermissionError: [Errno 1] Operation not permitted

Run log

RECORDED OUTPUT

Run failed: PermissionError: [Errno 1] Operation not permitted

Exception ignored in: <function MonitorClient.__del__ at 0x7bdcca2e2440>

Traceback (most recent call last):

```
File "/home/augustus/Documents/code/repos/pyrpl/pyrpl/redpitaya_client.py", line 94, in __del__
    self.close()
```

```
File "/home/augustus/Documents/code/repos/pyrpl/pyrpl/redpitaya_client.py", line 88, in close
    self.socket.send(b"c" + bytes(bytearray([0, 0, 0, 0, 0, 0, 0, 0])))
```

AttributeError: 'MonitorClient' object has no attribute 'socket'

Active parameters

RP commands: volts, bounded to 0–1. V_lambda: terminal volts for 2π . Gains: terminal V / RP V.

Parameter	Recorded value
label	"pd-visibility"
duration_s	12.0
report	"pdf"
pd_input	"in2"
visibility_frequency_hz	0.5
visibility_dark_voltage_v	null